

# Fetal Protocols

## Body

T2 Haste Coronals (mom uterus)

FOV 300  
Slice 5.0  
Dist. 20  
TR 1100  
TE 76  
Matrix 256 at 75%

T2 Haste Sagital (mom uterus)

FOV 300  
Slice 5.0  
Dist. 20  
TR 1100  
TE 76  
Matrix 256 at 75%

T2 Haste Sagital (mom cervix)

FOV 240  
Slice 5.0  
Dist 20  
TR 1100  
TE 76  
Matrix 256 at 75%

T2 Haste 78 (3 planes)

FOV 300  
Slice 3.0  
Dist. 0  
TR 1100  
TE 78  
Matrix 256 at 75%

Trufi Breath Hold (3 planes)

FOV 300  
Slice 4.0  
Dist. 0  
TR 4.05  
TE 1.65  
Matrix 256 at 100%

Trufi3d(cor uterus only)

FOV 400

Slice 2.00

Dist. 20

TR 3.99

TE 1.74

Matrix 256 x 80%

T1 Flash (sag)

FOV 300

Slice 4.0

Dist. 0

TR 202.0

TE 4.76

Matrix 256 x 65%

VIBE (cor uterus only)

FOV 300

Slice 2.0

Dist. 20

TR 3.52

TE 1.2

Matrix 256 x 75%

EPI 2d (3 planes)

FOV 334

Slice 3.0

Dist. 0

TR 5290

TE 50

Matrix 256 x 75